

Integration has never been so easy.

Now, A-dec® delivery systems with 300 Deluxe and 500 Deluxe Plus touchpads accommodate the KaVo® KL 703 LED electric micromotor. The streamlined design and intuitive user interface make the system easy to operate.

Enjoy an exceptional practice experience, with premier equipment that complements your professional expertise and reflects your commitment to superior patient care. Elevate your practice with A-dec and KaVo.



KaVo KL 703 Electric Motor: Smart. Efficient. Streamlined.



INTRA™ LUX KL 703 LED Electric Motor



LIGHTWEIGHT FEEL. COMPACT DESIGN. ERGONOMIC SHAPE.

Comfort and performance come together in the KaVo KL 703. Lightweight and compact, this sterilizable motor provides constant torque over the entire speed range for fatigue-free comfort and endless power. Suitable for any A-dec 200, A-dec 300 and A-dec 500 delivery system with a Deluxe or Deluxe Plus touchpad.

Features and Benefits

- Seamless electric motor integration with the choice between single and dual motors
- High-speed, low-speed, and torque-controlled rotary and reciprocating endo applications
- Speed range between 100 and 40,000 rpm with a 1:1 transmission ratio
- Brushless and sterilizable
- Up to 8 speed memory presents available in prep mode
- Up to 16 programmable memory settings (with A-dec 500 Deluxe Plus touchpad)
- ISO 6964 E-style coupling

Applications	General restorative, prophylaxis, endodontics	Length (mm)	34
A-dec Delivery System	A-dec 500B / A-dec 300 / A-dec 200	Weight (gm)	69
System	E-style	Sterilizable	Yes
Power (watts)	45	Micromotor/Type	Brushless motor
Speed (rpm)	100 to 40,000	Warranty	3 years
Torque (Ncm)	Max 3.0		

ELECTRIC MOTOR SYSTEM

Control Module + Motor + Tubing Kit			Motor + Tubing Kit	
Control Module			Motor	
	EMC-703 electric motor controller			
Motor			Tubing Kit	
	INTRA LUX KL 703 LED electric motor			
			Motor	
			Tubing Kit	
1st Position			2nd Position	